



Ultra-Stiff C+ Files.

Ideal for Initial Instrumentation and Tough Calcification.

Introducing stiff competition for extreme calcification – premium C+ Files from Maillefer. With their pyramid-shaped tip and unique taper for amazing stiffness and strength, these premium files make it easier to locate and access challenging canals. And because they are proven to provide up to a 300% gain in buckling force over standard K-Files, you'll be able to negotiate even the toughest calcification with confidence.

Strong, Superior Performance.

- Unique taper for ultra-stiff design
- Pyramid-shaped tip for easier penetration
- Highly polished for debris removal and smooth insertion into the canal
- Square section for resistance to distortion and increased safety
- Two available lengths with consistent cutting sensation across the whole range of instruments
- Calibration rings for easier depth determination

DENTSPLY
MAILLEFER

In Endodontics, Quality Has A Name.



MAILLEFER

DENTSPLY MAILLEFER ENDODONTICS

C+ Files

Significant Gain in Buckling Force.

- Easier location of canal orifices
- Easier access to apical third
- Better negotiation of calcifications
- Better tactile feel



Packages of 6 C+ Files

Size	21mm	25mm
06	670900	670904
08	670901	670905
10	670902	670906
15	670903	670907

C+ Files Taking Stiffness and Strength to the Next Level.

Get the premium file expertly designed to answer the challenge of highly calcified canals. Ultra-stiff C+ Files from Maillefer. Unwavering strength when you need it most.