


Gripping
Performance.

**The Senseus™ FlexoFile®
with a Premium, Soft-Grip Handle.**

The Senseus™ FlexoFile® instrument's unique, silicone handle is as intelligently designed as its flexible file. With a soft-to-the-touch, ergonomic design and improved tactile feedback, it's the next big thing in endodontic instrumentation.

senseus™
FlexoFile® Endodontic Instruments
The New Endo Sensation

DENSPY
MAILLEFER
In Endodontics, Quality Has A Name.
MAILLEFER



senseus™
FlexoFile® Endodontic Instruments
The New Endo Sensation

Large, Highly Ergonomic Handle Made From Soft-Texture Silicone

- Greater working comfort.
- Improved tactile feedback.
- Improved torque and force transmission.
- Excellent grip when used with gloves.
- X-ray visible calibration rings for accuracy and improved diagnostics.

Flexible, Stainless Steel File

- The twisted file provides flexible strength for preparing even the most curved and narrow canals.
- Safe-ended Batt® Tip.
- Preinserted silicone stops save you time and money.

Packages of 6 Senseus FlexoFile Instruments.

ITEM	SKU
21mm 06	670615
21mm 08	670616
21mm 10	670617
21 mm 15	670618
21mm 20	670619
21 mm 25	670620
21mm 30	670621
21mm 35	670622
21mm 40	670623
21mm 45	670624
21mm 15-40	670625
21mm 50	670626
21mm 55	670627
21mm 60	670628
21mm 70	670629
21mm 80	670630
21mm 45-80	670631
25mm 06	670632
25mm 08	670633
25mm 10	670634
25mm 15	670635
25mm 20	670636
25mm 25	670637
25mm 30	670638
25mm 35	670639
25mm 40	670640
25mm 15-40	670641
25mm 45	670642
25mm 50	670643
25mm 55	670644
25mm 60	670645
25mm 70	670646
25mm 80	670647
25mm 45-80	670648
31mm 06	670649
31mm 08	670650
31mm 10	670651
31mm 15	670652
31mm 20	670653
31mm 25	670654
31mm 30	670655
31mm 35	670656
31mm 40	670657
31mm 15-40	670658